

RENO-S90

More Great Sterilization

>> Specification

TOTAL CYCLE TIME

Non-Lumen Cycle About 28min.

Eco Cycle About 37min.

Advanced Cycle About 53min.

PROCESSING TEMPERATURE

Below 60°C

DIMENSIONS & WEIGHT

Dimensions 1,550(H)x700(W)x985(D)mm

Weight 390Kg

STERILIZATION CHAMBER

Configuration Rectangular

Total Volume 90.4L, 340(H)x380(W)x700(D)mm

Shelf Removable Two-tiered Shelf

Material Stainless Steel

DOOR

Door Type Single Door

Door Control Manual

STERILANT

Liquid Chemical Agent Hydrogen Peroxide

Concentration & Usage 8ml (50%), 1 Cycle / Cassette

CONNECTION

Electricity 230V~, Single Phase, 50/60Hz

or 200V Single Phase, 50/60Hz

Power Consumption 3 KVA

PROGRAM

Control PLC (Network)

Self-Test Mode Automatic

Printer Outputs Built-in Thermal Printer

Data Storage SD Memory Card / LAN / USB

DISPLAY

Screen Wide Touch Screen

Alarm Audible & Visible Alarm

MOBILITY

Mobile wheels or Cart Built-in 4 locking wheels

AIR INLET BACTERIA RETENTION FILTER

HEFA Filter 99.97% Efficiency at 0.3 micron (0.3×10^{-6})

REGULATORY APPROVALS

CE / ISO 13485 / ISO 9001

OPERATIONAL ENVIRONMENT

Temperature 5 ~ 40°C

Humidity 5°C~30°C / below 80%, 31°C~40°C / below 50%



Distributed by:

PT. Megah Alkesindo

Kompleks Ruko Rawa Bambu

Jl. Rawa Bambu Raya No. 15 G-H

Pasar Minggu Jakarta Selatan 12520 Indonesia

Phone : +62 21 788 44990 (Hunting)

Service : +62 21 2899 4111

Fax : +62 21 788 47144

Email : info@megahalkesindo.com