

RENO

Low Temperature Plasma Sterilizer

Expert in Infection Control




RENO-S30

“More reliable sterilization”

RENO-S30 Specification

Item	RENO-S30
External dimension	1,020(H) x 570(W) x 840(D)mm
Weight	190 Kg
Chamber	Dimension
	300(W) x 185(H) x 610(D)
Sterilant	Volume & Type
	34 L Rectanguler
Sterilant	Dose
	4ml (50%)
Sterilant	Cassette
	1 cycle/ cassette
Sterilant	Processing temperature
	Below 55°C
Sterilant	Total cycle Time
	About 45min.
Sterilant	Power Consumption
	2.5KVA
Sterilant	Electricity
	230V-Single Phase, 50/60Hz
Sterilant	Control
	Micro Processor
Sterilant	Display
	Touch Panel
Sterilant	Printer
	Thermal Printer
Sterilant	Mobility
	4 wheels



Color-coded touch screen



Sterilant residual-free system



Single-use cassette type sterilant

Distributed by:

PT. Megah Alkesindo

Kompleks Ruko Rawa Bambu

Jl. Rawa Bambu Raya No. 15 G-H

Pasar Minggu Jakarta Selatan 12520 Indonesia

Phone : +62 21 788 44990 (Hunting)

Service : +62 21 2899 4111

Fax : +62 21 788 47144

Email : info@megahalkesindo.com



RENOSEM

RENOSEM Co., Ltd

Rm-303-101, Techno-Park III,

36-1 Samjeong-dong, Ojeong-gu, Bucheon-shi

Gyeonggi-do, (#421-742) Korea