

SPESIFIKASI

NO	PRODUK	SPESIFIKASI	MERK	TYPE	NEGARA ASAL
1	Low Temperature Plasma Gas Sterilizer	MAIN CHARACTERISTICS • Colored wide touch screen control panel which displays cycle information • Control program by microprocessor • Spacious rectangular chamber (usable volume more than 90%) • Easy-to-clean chamber wall • Freely contactable chamber wall • Sterilant residual-free using vaporizing filter system • Built in Thermal Printer • Plasma generated by dielectric barrier discharge & corona plasma • Provided with data backup port (using Memory Card) • Sterilization process with two ways reactions (Radical Plasma & Abatement Plasma)	RENOSEM	RENO-S30	Korea
		GENERAL INFORMATION • Sterilizing temperature : Below 55°C • Sterilizing pressure: 500mTorr - 20Torr			
		CHAMBER • Chamber volume 34 Liter • Dimension (mm) : Approx. 300 x 185 x 610 • Made of high quality stainless steel			
		COMPATIBLE ITEMS • Metal materials (Stainless Steel, Aluminum, Titanium alloys etc.) • Plastic products (Polyethylene, Polycarbonate, PVC, PP, Viton etc.) • Glass materials • Rubber materials (Silicon) • Ceramics material • Stainless steel lumen: Diameter >1ø mm and up to 1.000 mm in length • Flexible lumen: Diameter >1ø mm and up to 12.000mm in length • Single through lumen: Diameter > 2ø mm and up to 1.500mm in length			
		STERILANT • Hidrogen Peroxide (50%) • 1 cassette/cycle			
		ELECTRICITY • Rated voltages: 200V/220V/230V; 50/60 Hz • Power consumption: 2,5 kVA • Max current 12A			