


RENO

Low Temperature Plasma Sterilizer



Reno-D50 Specification

Item	Reno-D50
External dimension	1,310(H) x 538(W) x 789(D)
Weight	290Kg
Chamber	Dimension: 185(H) x 300(W) x 550(D) Volume & Type: 60L, 2 x Chamber(30L), Rectangular
Sterilant	Dose: 4ml x 2chambers (50%) Cassette: 1 Cycle/2cassette
Processing temperature	Below 55°C
Total cycle Time	About 45min
Power Consumption	2.5KVA
Electricity	230V-Single Phase, 50/60Hz
Control	Micro Processor
Display	Touch Panel
Printer	Thermal Printer
Mobility	4 wheels



Color-coded touch screen



Sterilant residual-free system



Single-use cassette type sterilant

Distributed by:

PT. Megah Alkesindo

Kompleks Ruko Rawa Bambu
Jl. Rawa Bambu Raya No. 15 G-H
Pasar Minggu Jakarta Selatan 12520 Indonesia
Phone : +62 21 788 44990 (Hunting)
Service : +62 21 2899 4111
Fax : +62 21 788 47144
Email : info@megahalkesindo.com



RENOSEM Co., Ltd
Rm-303-101, Techno-Park III,
36-1 Samjeong-dong, Ojeong-gu, Bucheon-shi
Gyeonggi-do, (#421-742) Korea