


ENT Ear & Nasal Irrigator CMEI-100

This is a portable ENT treatment for ear and nose. Water Temperatur (1-45°C) and Operation time (5-90 seconds) can be adjusted by user preference. Ear and nose tip can be changed for Hygienic treatment as well. As adophth enough 1,5 liter tank. It can be used as portable irrigator or attach on ENT unit.



Irrigator For Ear and Noise CMEI-100

Specification	
Power source	AC 220V, 50/60Hz
Power Consumption	310VA
Temperature	1°~45°
Use time	5~90(SEC)
Dimension	160mm(W)x330mm(D)x331mm(H)
Weight	10Kg
Standart unit	Main body
	Storage
	Funnel
	Irrigator hand gun
	Ear irrigator
	Nasal Irrigator 90° Tip
	Irrigator bottle (1500CC)

