

SPESIFIKASI

NO	PRODUK	SPESIFIKASI	MERK	TYPE	NEGARA ASAL
1	CHAMMED Irrigator For Ear and Noise With Arm Set - Chammed CMEI-100	<p>Features</p> <p>This is a portable Ent treatment for Ear and Nose. Water Temperature(1~45°C) and operation time (5~90 seconds) can be adgusted by users preference. Ear and nose tip can be changed for Hygienic treatement as well. As adopt enough 1.5 Liter water tank, it can be used as a portable irrigator or attach on ENT unit. For patient's satety, safety systems are applied on water tank cover, heater, temperature sensor and water level connector.</p> <p>Specification</p> <p>Power source : AC 220V, 50/60Hz</p> <p>Power Consumption : 310VA</p> <p>Temperature : 1°~45°</p> <p>Use time : 5~90(SEC)</p> <p>Dimension : 160mm(W)x330mm(D)x331mm(H)</p> <p>Weight : 10Kg</p> <p>Standart unit : Main body</p> <p>Storage</p> <p>Funnel</p> <p>Irrigator hand gun</p> <p>Ear irrigator</p> <p>Nasal Irrigator 90° Tip</p> <p>Irrigator bottle (1500CC)</p>	Chammed	CMEI-100	Korea