ENT Treatment CU-5000 SPESIFIKASI

| NO | PRODUK | | MERK | TYPE | NEGARA ASAL |
|----|------------------|--|---------|---------|----------------|
| 1 | CHAMMED ENT | ENT Treatment CU-5000 | CHAMMED | CU-5000 | KOREA |
| | Treatment Unit - | - It's full system for treatment and operating in Otolaryngology. | | | |
| | Chammed CU-5000 | - Efficient 2-suctions | | | |
| | | - Convenient total control switches can control each functions easily | | | |
| | | - Safety auto-off switch for Heater, either manual switch or timer control switch(30sec) | | | |
| | | - Glass suction bottles are enough size in usage and good in cleaning | | | |
| | | - Artificial marble surface enables dust proof, easy to clean | | | |
| | | - Drawer typed built in Endoscope storage & UV sterilizer optimizes usability | | | |
| | | Various optional accessories enable enlarge user's choice | | | |
| | | - Personal computer rack | | | |
| | | - Key board supporter with swiveling 2-folded arm | | | |
| | | - Marble top can use desk, chart or multi-function | | | |
| | | - Front body 105mm stainless cans, which are convenience for used instrument container | | | |
| | | - Suction bottle supports are drawer type which are enable to clean easily | | | |
| | | - Two-suction system enable user choose either one | | | |
| | | - Maximum 3 post pole to add up accessories conveniently | | | |
| | | - Six 83mm stainless can give enough space to use | | | |
| | | - Casters can move device easily to other place | | | |
| | | - Device surface easy to clean and keep disinfection | | | |
| | | - Color variety: Gray or white | | | |
| | | - Main Instrument tray on top of device has dust proof cover and can heat for instrument | | | |
| | | - Suction bottle window enable user can check bottle liquid level easily | | | |
| | | - Medicine bottles with folded arm supporter are convenient to use | | | |
| | | - Suction device is automated micro-switch movement | | | |
| | | - Sprays are working automatically when pick up and stop working return to the holder | | | |
| | | - Instrument tray cover can protect instrument from being spoiled by dust and should semi- | | | |
| | | transparency to check inside also | | l | |
| | | Instrument cover should be made of acryl, which is not broken easily | | | |
| | | SPECIFICATIONS ; | | | |
| | | - Power source: AC220V, 50/60Hz | | | |
| | | - Power consumption: 1200VA(Main body) | | | |
| | | - Suction motor: 100ℓ/min, 680mmHg | | | |
| | | - Compressor motor 2.3Kgf/ait±0.3Kgf/ait | | | |
| | | - Dimensions: 1800(w) * 709(d) * 800(h), Post pole max. height: 1946mm(diff. 10%) | | | |
| | | - Weight : Approx 178 kg ± 10 % | | | |
| | | Main body & standard accessories | | | |
| | | - Illumination light (halogen reflector lamp, R80, Flood) | | | |
| | | - Medicine bottles (Transparent, brown,blue) with tray | | | |
| | | - Spray: 4(Straight 3, Curved 1) | | | |
| | | - Suction: 2 ea | | | |
| | | - Nasal Suction Tip (# 1-4) | | | |
| | | - Suction bottle 2500cc / 1000cc / 500cc, each 2ea | | | |
| | | - 100mm : 2ea / 83mm : 9ea / 53mm: 3ea Can | | | |
| | | - Instrument Tray with cover, (with Pre-heating) | | | |
| | | - Anti Fog (quick mirror warmer) with auto cover open / close | | | |
| | | - Auto politzer (ventilation) for Eustachian tube with tips (10ea) | | | |
| | | - UV Sterizer for scopes (drawer type) | | | |
| | | - Side can set | | | |
| | | - Chair control switch (Up / Down / Rear / Front / Reset) | | | |