INSPITAL



Operating Table OT80.20

NSP!TAI

Feature



Functional Design & Wide Range of Accessory **Options**

- · Modular design taking shape by different surgical operations with high stability
- Solutions ensuring to reach the position required by any kind of surgical operations with its wide range of accessory options



Sliding

· Image reception ability in a wide area without disturbing the sterile environment by its sliding feature within 300 mm on horizontal plane



Optimal Clearance In Surgical Area

· Maximum surgical area visibility with its kidney elevator, integrated on the table and adjustable up to1500(70 mm)



Easy Approach for C-Arm

- · Complete access to all table areas and its X-ray permeable table top allow high image quality
- Suitable base structure for easy access and approach of the C-arm scopy and surgical microscopes



T-Base Specifications

- · T-Base enables table optimal approach for the surgical team and equipment
- · ABS Material with break and corrosion resistance
- · Subbase is made of ABS plastic material. Wall thickness is 5 mm.



Optimum Mattress Performance

· Inspital operating table mattresses are designed and manufactured according to the Pressure Mapping test results.



Different Mattress Options

- · Compact material table top with high X-ray permeability that allows high-quality image reception
- · Visco or polyurethane mattress options that ensure the patient's comfort by taking its body shape, thus preventing decubitus prophylaxis



Accessory Positioning

· Stainless steel rail attachment system integrated on both sides of the table, enabling accessory mounting without removing existing accessories

TECHNICAL FEATURES

· Sliding Movement : 300 mm · Lateral Tilt Movement : +20°/-20° · Back Rest Movement : +80°/-45° Trendelenbura $: +30^{\circ}$ Reverse Trendelenburg : -30° : +30°/-80° Leg Rest Up / Down : +30°/-90° · Head Rest Movement

• Flex / Reflex Position : + 225° / - 100° · Width x Length · Height (min) Height (max) Mattress Thickness

 Castors Diameter · Max Patient Weight Kidney Surgery Position

· Stainless Steel T-base · Visco mattress

: 560x2200 mm

: 595 mm : 1035 mm 70 mm Ø 100 mm

Back Rest

360 kg / 800 lbs. 150° (optional), OT80.20 Optional, OT80.20S

Optional, OT80.20V

 Voltage : 220-240V AC 50/60Hz

: 24 VDC Working Voltage

· Fuse Type : 2 x 1,25A Glass Fuse Battery Duration : 30 Full Movement


 Battery type : 2 x 12 V · Battery Charge Time · 9 Hour

: Class I (MDD 93/42 EEC) Classification

 Standards : TS EN60601-2-46

ve MDD/93/42/EEC


Electro-hydraulic



Manual



Leg Rests, Up / Down



Head Rest



Operating Table OT80.20

INSPITAL

ACCESORIS OPERATING TABLE



OT60.90 Orthopedic Traction Set



OT60.82 Cross Bar Adapter



OT60.80 SET Neurosurgery Head Rest with Pin &



OT60.76 Arthroscopy



OT60.77 Leg Positioning



OT60.01 Arm Support,



OT60.03 Anaesthesia



OT60.07 Body Fixing Strap



OT60.05 IV Stand, 2 hooks



OT60.06 IV Stand, 4 hooks



OT60.04 Anaeshtsia



OT60.71 Patient Heating



OT60.08 Geopel Leg Rest



OT60.11 Irrigation Bowl,



OT60.12 Irrigat. Bowl With



OT60.14 Waste Bin



OT60.13 Irriga. Bowl With



OT60.17 Clamp, Rotatable

& Adjustable Jaw



OT60.18 Clamp, Rotatable



OT60.15 Clamp, Easy Fix



OT60.16 Clamp



OT60.09 Shoulder Support



OT60.10 Lateral Support



OT60.19 Foot Support



OT60.20 Proctology Set



OT60.21 Ophtalmologic



OT60.22 Weinberger Hand



OT60.23 Leg Strap



OT60.24 Meniscus Posi-



OT60.72 Knee Arthroscopy



OT60.73 Carbonfiber



OT60.25 Sliding Handles



OT60.26 Arm Surgery



OT60.29 Support for Hand



OT60.30 X-Ray Cassette Support



OT60.31 X-Ray Cassette Tray



OT60.33 Orthopedic Traction Set Trolley



OT60.32 Accessory Trolley



BS10.20 Step, Double


OT60.54 Arm Rest, Gel




BS10.10 Step, Single



OT60.61 Wilson Frame Set




OT60.62 Gel Mattress, Adult






OT60.50 Head Rest, Gel, Adult OT60.40 Head Rest, PU, Adult



OT60.59 Heel Support, Gel



Gel. Adult OT60.53 for Pediatric



OT60.58 Foot Wrist Support, Gel




OT60.4 Horseshoe Head Rest, PU, Adult



OT60.60 Lateral Positioner, Gel



OT60.57 Prone Head Support, Gel




OT60.55 Chest Roll, Gel,

Adult



OT60.42 Chest Roll, PU



Operating Table OT80.20

INSPITAL



Left Lateral Tilt Position



Struma Surgery Position



Reflex Position



Flex Position



Chest (Breast) Surgery Position Chest Surgery Position





Gynecologic Surgery Position



Spinal Surgery Position



Meniscus Surgery Position



Ophthalmology Position



Urologic Surgery Position



Orth.Fracture Surgery Position



Leg Positioning Boot



Leg Positioning Boot







- ♦ Made of 304 quality stainless steel
- ♦ Leg strainer function by mechanical lever
- ◆ Traction boots with up & down movement and rotation capability around its own axis



Neurosurgery Head Rest with Pin & Horseshoe Head Rest Designed for use in neurosurgery

- ◆ Extremely durable and rigid construction from aluminum alloy
- ◆ 3 pcs fixing system with pins & completely adjust. skull clamps
- ◆ Adjustment mechanism from 4 different points
- Easily mountable horse shoe head rest accessory A.C: OT60.80 SET
- \blacklozenge Can be used with skull fixing pins or horseshoe head rest accessories





INSPITAL INSPITAL Medical Technology GmbH

Karaoğlan Mahallesi Küme Evleri No: 745 Gölbaşı / ANKARA, TÜRKİYE Tel: +90 312 619 02 22 Fax: +90 312 619 02 25

REG. NO: KEMENKES RI AKL.....

PT MEGAH ALKESINDO

Jl. Rawa Bambu Raya No. 15 F-H

Pasar Minggu Jakarta Selatan 12520 Indonesia Phone: +62 21 788 77990 (Hunting) Fax: +62 21 788 47144

 ${\sf Email:info@megahalkesindo.com}\ \ {\sf Website:www.megahalkesindo.com}$